



New Era Pump Systems Inc.

www.SyringePump.com

NE-8060 Syringe Pump

Model NE-8060 High Pressure / High Speed Syringe Pump

Features

- All Programming features of the NE-1000 syringe pump.
- Syringes size up to 140 mL. 200 mL partially filled (About 160 mL).
- Pumping force over 200 lbs.

Adjustable force limit: Force adjustment knob reduces the maximum force to reduce damage to syringes.

Limit Switches

- Two adjustable limit switches: One sets the infusion limit, one sets the withdrawal limit. To adjust the limit switches, loosen the white thumbscrew, then slide the collar to the required position. Then re-tighten the thumbscrew.
- When the pusher block reaches a limit switch, the pump will either: 1) Trigger a Program Event function; 2) Continue with the next Program Phase; 3) Stop the pump.

If a Program Event is set (Event Function), the Event trap will be triggered and the Pumping Program will continue execution with the Program Phase set in the Event Function.

Otherwise, the Pumping Program will start the next Program Phase unless the next Program Phase pumps in the same direction. This is a safety feature to prevent attempting to pump past the limit switch. The pump will instead stop.

Limit switches can be used as a program volume target. Limit switches will override the volume target setting.

Additional Notes

- Does not have stall detection.
- Pusher block has a solid nut block. It is not releasable.
- To position the pusher block, use the purge function: Set the pumping direction with the direction key '→←', then press and hold the 'Start/Stop' key. The display will indicate **PURG** and the pusher block will move at top speed in the set direction. Release the 'Start/Stop' key to stop the pump.

Maintenance

- Periodic lubrication is required for proper operation. **Failure of the nut block can occur if not properly lubricated.**
- Lead-screw: Grease
- Guide Rods: Oil
- Remove accumulation of dirt or debris

WARNING

Use extreme caution. The NE-8060 can deliver enough force to break syringes or any objects that get caught in the mechanism. Broken syringes and associated plumbing can create dangerous projectiles that can cause bodily harm. Keep fingers and loose clothing away from mechanism.

Syringe Manufacturer (all names ™)	Syringe (mL)	Inside Diameter (mm)	Minimum Rate (µL/hr)	Maximum Rate (mL/min)
B-D	1	4.699	6.803	15.62
	3	8.585	22.71	52.15
	5	11.99	44.3	101.7
	10	14.43	64.16	147.3
	20	19.05	111.9	256.7
	30	21.59	143.7	329.8
	60	26.59	217.9	500.2
HSW Norm-Ject	1	4.69	6.777	15.56
	3	9.65	28.7	65.89
	5	12.45	47.76	109.6
	10	15.9	77.89	178.8
	20	20.05	123.9	284.4
	30	22.9	161.6	371
	50	29.2	262.7	603.3
Monoject	1	5.74	10.16	23.31
	3	8.941	24.63	56.56
	6	12.7	49.7	114.1
	12	15.72	76.14	174.8
	20	20.12	124.8	286.4
	35	23.52	170.5	391.4
	60	26.64	218.7	502.1
	140	38	444.9	1021
Terumo	1	4.7	6.806	15.63
	3	8.95	24.68	56.67
	5	13	52.07	119.5
	10	15.8	76.92	176.6
	20	20.15	125.1	287.2
	30	23.1	164.4	377.5
	60	29.7	271.8	624.1
Poulten & Graf (Glass)	1	6.7	13.84	31.76
	2	8.91	24.46	56.17
	3	9.06	25.29	58.08
	5	11.75	42.54	97.69
	10	14.67	66.31	152.2
	20	19.62	118.6	272.3
	30	22.69	158.7	364.2
	50	26.96	224	514.3
Steel Syringes	1	9.538	28.03	64.37
	3	9.538	28.03	64.37
	5	12.7	49.7	114.1
	8	9.538	28.03	64.37
	20	19.13	112.8	258.9
	50	28.6	252.1	578.7
	100	34.93	376	863.3
	200	44.75	617	1416

SGE Syringe (mL)	Inside Diameter (mm)	Maximum Rate (mL/hr)	Minimum Rate (µL/hr)	SGE Syringe (mL)	Inside Diameter (mm)	Maximum Rate (mL/min)	Minimum Rate (µL/hr)
0.25	2.303	225.1	1.635	10	14.57	150.2	65.41
0.5	3.257	450.3	3.269	25	23.03	375.2	163.5
1	4.606	900.6	6.537	50	27.5	535.1	233
2.5	7.284	2252	16.35	100	34.99	866.2	377.2
5	10.3	4504	32.69				

Specifications

<u>Model</u>	<u>Style</u>	<u>Stall Detection</u>	<u>Number of Syringes</u>	<u>Maximum Syringe Size</u>
NE-8060	Stand-Alone	No	1	140 mL 200 mL (Partially filled 160 mL)

Mechanical

Drive block type: Solid (Must use Purge function to move pusher block).

Motor type: Step motor
Motor steps per revolution: 200
Motor to drive screw ratio: 1/1
Drive screw pitch: 8 revolutions/”

Dimensions: 11 1/4” x 6 1/8” x 6 3/8” LxWxH
(28.575 cm x 15.5575 cm x 16.1925 cm)

Weight: 9.5 lbs. (4.31 kg)

Electrical

Power supply type: External, country and power source specific
Power supply output rating: 24V DC @ 2.5A
Power connector: 2.1 mm, center positive, DC

Operational

Accuracy: Within 1% error
Reproducibility: Within 0.1% error
Maximum force: 200 lbs.
Micro-stepping: 1/8 to 1/1 depending on motor speed
Advance per step: 1.984375 μ m to 15.875 μ m depending on motor speed
Maximum speed: 90.12566735 cm/min
Minimum speed: 0.039242 cm/hr
Maximum pumping rate: 1021 mL/min with a 140 mL syringe
Minimum pumping rate: 6.803 μ L/hr with a B-D 1 mL syringe
Syringe inside diameter range: 0.100 to 50.00 mm
Number of Program Phases: 41
RS-232 pump network: 100 pumps maximum
RS-232 selectable baud rates: 300, 1200, 2400, 9600, 19200

