



Disposable Syringes & Kits

New Era supplies an assortment of disposable plastic syringes between 1mL to 140mL from Monoject, NormJect, Terumo & Excel. You can receive discounts in box quantities.

Description (mL=CC)	Part #	Price Each	Box Part #	Box Price
140 mL Luer Lock Non-Inert	SYR-MJ140LL	\$6.00	SYR-MJ140LLBX	\$119.00 (20 per Box)
60 mL Luer Lock Non-Inert	SYR-TER60LL	\$1.25	SYR-TER60LLBX	\$18.75 (25 per Box)
50 mL Luer Lock Inert	SYR-NJ50LL	\$1.75	SYR-NJ50LLBX	\$37.50 (30 per Box)
30 mL Luer Lock Inert	SYR-NJ30LL	\$1.25	SYR-NJ30LLBX	\$40.00 (50 per Box)
20 mL Luer Lock Inert	SYR-NJ20LL	\$0.75	SYR-NJ20LLBX	\$55.00 (100 per Box)
10 mL Luer Lock Inert	SYR-NJ10LL	\$0.50	SYR-NJ10LLBX	\$40.00 (100 per Box)
5 mL Luer Lock Inert	SYR-NJ5LL	\$0.40	SYR-NJ5LLBX	\$30.00 (100 per Box)
3 mL Luer Lock Inert	SYR-NJ3LL	\$0.35	SYR-NJ3LLBX	\$25.00 (100 per Box)
1 mL Luer-Lock Non-Inert	SYR-EX1LL	\$0.40	SYR-EX1LLBX	\$30.00 (100 per Box)
1 mL Tuberculin Inert	SYR-TB1	\$0.25	SYR-EX1LLBX	\$20.00 (100 per Box)
Large Syringe kit	P-SYRKIT-LG	\$15	N/A	N/A
Small Syringe kit	P-SYRKIT-SM	\$11.50	N/A	N/A
Dual Pump Plumbing Kit	P-DKIT	\$20.00	N/A	N/A



Large Syringe Kit
Part ID: P-SYRKIT-LG

- (2) 20mL Luer Lock Inert
- (2) 30mL Luer Lock Inert
- (2) 60mL Terumo Luer Lock Non-Inert
- (6) Female Luer Barbs
- (3) Male Luer Barbs
- (20 ft) Tubing .125" ID; .188" OD, PVC



Small Syringe Kit
Part ID: P-SYRKIT-SM

- (2) 1mL Luer Lock Inert
- (2) 3mL Luer Lock Inert
- (2) 5mL Luer Lock Inert
- (2) 10mL Luer Lock Inert
- (6) Female Luer Barbs
- (3) Male Luer Barbs
- (20 ft) Tubing .125" ID; .188" OD, PVC



Dual Pump Plumbing Kit
Part ID: P-DKIT

- For use with our continuous infusion system or dual infusion system (Dual-NE-1000). Provides a dual check valve for two syringes, withdraw will pull fluid from a reservoir and infusion will output into a separate line
- The kit is pre-assembled, all you need to do is plug in the syringes
- Each dual check valve uses 6" of 1/8" tubing to connect to a Y connector, one for reservoir and one for output; each Y-connector has 5' of tubing to make connections to your

application. Two 60cc Terumo syringes are included

- Highly corrosive chemicals are not compatible with this kit. Contact us for a solution that is compatible.